


# Enlight

5th Generation of Homogeneous CLIA




Let's enlighten diagnostics  
together with "Enlight"!

## Working Principle of CRET (Chemiluminescence Resonance Energy Transfer)

The Chemiluminescence Resonance Energy Transfer (CRET) method is based on the energy transfer between a donor (antigen) molecule and an acceptor (antibody) molecule. This energy transfer enhances the chemiluminescent signal and is employed in advanced diagnostic applications to achieve high sensitivity, precision, and efficiency.

### Chemiluminescence Resonance Energy Transfer (CRET) Principle



### STEPS OF CRET

- 1. Chemical Reaction:**  
A chemical reaction excites a donor (antigen) molecule, giving it energy.
- 2. Energy Transfer:**  
The donor (antigen) molecule passes this energy to a nearby acceptor (antibody) molecule. This happens through a process called resonance energy transfer (like a wireless energy transfer between two molecules).
- 3. Light Emission:**  
The acceptor (antibody) molecule uses this energy to emit a bright light signal, which is stronger and easier to detect than the donor's original light.
- 4. Signal Detection:**  
The emitted light is detected and measured, and its intensity tells how much of the target analyte (e.g., a hormone or protein) is present.

### Advantages of CRET

- High Sensitivity : Amplified signals allow detection of low-concentration analytes.
- Fast Response : Minimal reaction steps make CRET-based assays rapid.
- Stable Signals : The energy transfer mechanism ensures a more stable and reliable output.
- Compact Design : CRET enables the development of portable diagnostic systems, suitable for POCT (Point of Care Testing).

## COMPARISON of CRET, CLIA, and ECLIA

FEATURE	CRET	CLIA	ECLIA
<b>Principle</b>	Resonance Energy Transfer (RET) between antigen and antibody molecules, enhancing the chemiluminescent signal.	Enzyme-substrate reactions producing light, proportional to the an alyte concentration.	Electrochemical excitation of luminescent molecules, generating light at an electrode
<b>Sensitivity</b>	Highest sensitivity due to amplified signal via RET.	High sensitivity; suitable for detecting low concentration an alytes.	Very high sensitivity; ideal for precise diagnostics.
<b>Speed of Analysis</b>	Fastest, with fewer reaction steps.	Moderate; enzymatic reactions may require washing and incubation steps.	Fast; electrochemical reactions are efficient and controlled.
<b>Signal Stability</b>	Highly stable due to donor-acceptor specificity.	Variable; depends on enzyme and substrate stability.	Extremely stable; electrochemical control ensures reproducibility.
<b>Cost (Setup)</b>	Moderate; advanced RET materials increase initial cost.	Affordable; uses well established reagents and methods.	Expensive; requires specialized electrodes and equipment.
<b>Cost (Operational)</b>	Cost-effective; minimal reagent usage and faster testing.	Moderate; reagent consumption may increase costs.	High; involves costly reagents and maintenance of electrodes.
<b>Portability</b>	Highly portable; compact design ideal for POCT.	Moderate portability; typically requires medium sized systems.	Low; requires large, laboratory-based systems.
<b>Applications</b>	POCT, rapid testing in clinics and rural areas.	General diagnostics in centralized labs with high throughput.	Specialized diagnostics for cancer markers, cardiac panels, and critical care
<b>Throughput</b>	Medium to high; optimized for smaller setups.	High; suitable for bulk testing in large labs.	High; designed for fully automated, high throughput systems.
<b>Maintenance</b>	Low; fewer mechanical components.	Moderate; periodic enzyme and reagent management required.	High; regular electrode maintenance and calibration needed.
<b>Example Machines</b>	Enlight DC80 Dry Chemiluminescence Analyzer	Abbott Architect i1000SR, Mindray CL900i	Roche Cobas e411, Siemens ADVIA Centaur

# Chemiluminescence Resonance Energy Transfer (CRET) **VS** Traditional Chemiluminescence Immunoassay (CLIA)

FEATURE	CRET METHOD	TRADITIONAL CLIA
Working Principle	<ul style="list-style-type: none"> <li>▶ Based on energy transfer between a donor molecule and an acceptor molecule.</li> <li>▶ The energy transfer occurs when the donor molecule is excited and emits energy, which is then transferred to the acceptor molecule.</li> <li>▶ The resulting chemiluminescent signal is enhanced by the resonance energy transfer, improving sensitivity and efficiency.</li> </ul>	<ul style="list-style-type: none"> <li>▶ Involves the use of a chemical reaction where a chemiluminescent substrate emits light as a by-product.</li> <li>▶ The light emission is directly proportional to the analyte concentration.</li> <li>▶ Typically uses enzymatic or electrochemical reactions to generate the luminescent signal</li> </ul>
Sensitivity and Accuracy	<ul style="list-style-type: none"> <li>▶ <b>Higher sensitivity</b> due to the enhanced signal from resonance energy transfer.</li> <li>▶ More accurate at detecting low concentration analytes because of its ability to amplify the luminescent signal.</li> </ul>	<ul style="list-style-type: none"> <li>▶ Highly sensitive, but may not reach the amplification levels of CRET.</li> <li>▶ Accuracy is dependent on enzyme efficiency and substrate purity.</li> </ul>
Speed of Analysis	<ul style="list-style-type: none"> <li>▶ Faster due to the direct and efficient energy transfer mechanism.</li> <li>▶ Requires fewer steps in signal generation, reducing assay time.</li> </ul>	<ul style="list-style-type: none"> <li>▶ Relatively slower, as it depends on enzymatic reactions or secondary processes for signal generation.</li> <li>▶ Some protocols involve washing and incubation steps, increasing total assay time.</li> </ul>
Stability and Reliability	<ul style="list-style-type: none"> <li>▶ Chemiluminescent signal is more stable due to the specific donor-acceptor pairing.</li> <li>▶ Less prone to interference from external factors like temperature or pH.</li> </ul>	<ul style="list-style-type: none"> <li>▶ Signal stability can vary based on the enzyme and substrate used.</li> <li>▶ More prone to variability due to external factors.</li> </ul>





# Chemiluminescence Resonance Energy Transfer (CRET) **VS** Traditional Chemiluminescence Immunoassay (CLIA)

FEATURE	CRET METHOD	TRADITIONAL CLIA
Cost and Complexity	<ul style="list-style-type: none"> <li>▶ Initial setup may be costlier due to advanced materials and specific donor-acceptor pairs.</li> <li>▶ Cost-effective in the long run due to lower reagent consumption and faster analysis.</li> </ul>	<ul style="list-style-type: none"> <li>▶ More affordable initially, as it uses established chemistries and reagents.</li> <li>▶ May require higher reagent volumes, increasing operational costs over time</li> </ul>
Compactness and Portability	<ul style="list-style-type: none"> <li>▶ Designed for <b>compact and portable systems</b>, making it ideal for POCT (Point of Care Testing).</li> <li>▶ Simplifies the hardware requirements for light detection.</li> </ul>	<ul style="list-style-type: none"> <li>▶ Generally requires more sophisticated laboratory equipment, making it less portable.</li> </ul>
Applications	<ul style="list-style-type: none"> <li>▶ Suitable for rapid testing and decentralized diagnostics, such as POCT devices in clinics or rural areas.</li> <li>▶ Ideal for assays where speed, sensitivity, and compactness are critical.</li> </ul>	<ul style="list-style-type: none"> <li>▶ Widely used in centralized labs for high throughput testing.</li> <li>▶ Preferred for established protocols and a wide range of assays.</li> </ul>



*"Fast, flexible, cost-effective and highly sensitive Enlight is the perfect CLIA Instrument for labs of any size."*

**Encore Biomedicals Pvt. Ltd.**

 R-548, 3rd floor, TTC Industrial Area, Rabale, Navi Mumbai - 400701.  +91 91368 84688  [encorebiomedicals@gmail.com](mailto:encorebiomedicals@gmail.com)  [www.encorebiomedicals.com](http://www.encorebiomedicals.com)